
Therapy guide Endomed 581ID & Endomed 582ID

Therapiewijzer
Therapy manual
Therapiefibel
Guide de thérapie
Guia de terapia

Copyright:

**Enraf
Nonius**

Enraf-Nonius B.V.

P.O. Box 12080

3004 GB ROTTERDAM

The Netherlands

Tel.: (0)10 - 20 30 600

Fax: (0)10 - 20 30 699

Article number 1495.786-41

December 2005

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A. PAIN RELIEF	
I. SUPERFICIAL PAIN (DIA.)	23
II. NEUROGENIC PAIN (DIA.)	23
III. ACUTE PHASE (MF)	23
IV. SUBACUTE PHASE (MF)	23
V. CHRONIC PHASE (MF)	24
VI. ACUTE PHASE (TENS)	24
VII. SUBACUTE PHASE (TENS)	24
VIII. CHRONIC PHASE (TENS)	24
Points IX - XII only apply to the Endomed 582 ID	
IX. DEEP PAIN: ACUTE PHASE (MF)	25
X. DEEP PAIN: SUBACUTE PHASE (TENS)	25
XI. DEEP PAIN: SUBACUTE PHASE (MF)	25
XII. MYALGIA (MF)	25
B. CIRCULATION IMPROVEMENT	
I. SUPERFICIAL CIRCULATION IMPROVEMENT (DIA.)	26
II. CIRCULATION IMPROVEMENT (TENS)	26
III. CIRCULATION IMPROVEMENT (MF)	26
Point IV only applies to the Endomed 582 ID	
IV. INTERMITTENT CLAUDICATION (TENS)	26
C. MUSCLE RELAXATION	
I. SINGULAR PROGRAMME (TENS)	27
II. SEQUENTIAL PROGRAMME (MF)	27
D. NERVE LESIONS	
I. MUSCLE STIMULATION (FARADISM)	28
II. MUSCLE STIMULATION (TENS)	28
III. MUSCLE STIMULATION: PHASIC FIBRES (INITIAL STAGE)	28
IV. MUSCLE STIMULATION: PHASIC FIBRES (FINAL STAGE)	29
V. MUSCLE STIMULATION: TONIC FIBRES (INITIAL STAGE)	29
VI. MUSCLE STIMULATION: TONIC FIBRES (FINAL STAGE)	29
E. MUSCLE TRAINING	
I. NEUTRAL	30
II. SLOW-TWITCH FIBRES	30
III. FAST-TWITCH FIBRES	30
IV. TOTAL TRAINING	31
F. MUSCLE STRENGTHENING	
I. SINGULAR PROGRAMME (TENS)	32
II. SEQUENTIAL PROGRAMME	32
III. SEQUENTIAL PROGRAMME (PHASIC FIBRES)	32
IV. SEQUENTIAL PROGRAMME (TONIC FIBRES)	33
V. SINGULAR PROGRAMME (MF)	33
VI. PROGRAMME BASED ON 'RUSSIAN STIMULATION'	33
Chapters G and H only apply to the Endomed 582 ID	
G. INFLAMMATION	
I. BURSITIS (MF)	34
II. CAPSULITIS (MF)	34
III. EPICONDYLITIS (TENS)	34
IV. TENDINITIS (MF)	34
H. URINARY INCONTINENCE	
I. STRESS INCONTINENCE	35
II. URGE INCONTINENCE	35
III. MIXED INCONTINENCE	35

I. SUPERFICIAL PAIN (DIA.)

1. General: Programme for the treatment of pain related disorders by using the milder current forms of diadynamic current as described by dr. Bernard.
2. Indication: Superficial, sharp and/or dull pain.
3. Treatment time: 8 minutes.
4. Remarks: Diadynamic current is a Direct Current. At higher current densities DC may cause tissue damage. Intensity should not exceed 0.2mA/cm² of effective treatment area. Preferably use thick pads/sponges (min. 1 cm).
5. Programme: 2 min. DF, 3 min. LP, automatic polarity reverse
3 min. LP
6. Sensation: Tickling, prickling followed by alternately light to moderate vibrations.
7. Intensity level: Clearly to strongly perceptible, preferably just below motoric threshold.

II. NEUROGENIC (DIA.)

1. General: Programme for the treatment of neurogenic pain by using the CP form of diadynamic current (CP = courant modulé en Courtes Périodes) and automatic polarity reverse.
2. Indication: Neuralgia, causalgia, post-herpetic neuralgia, phantom limb pain (besides pain relief, local circulation is increased).
3. Treatment time: 6 minutes.
4. Remarks: Diadynamic current is a Direct Current. At higher current densities DC may cause tissue damage. Intensity should not exceed 0.2 mA/cm² of effective treatment area. Preferably use thick pads/sponges (min. 1 cm).
5. Programme: 3 min. CP, automatic polarity reverse, 3 min. CP
6. Sensation: Alternately light and rough vibrations.
7. Intensity level: Clearly perceptible, not painful. Preferably just below motoric threshold.

III. ACUTE PHASE (MF)

1. General: Programme for Medium Frequency current (interferential current, bipolar mode) using varying parameters in order to prevent accommodation.
2. Indication: Acute pain in absence of muscular hypertonia or as initial treatment for patients who are extremely sensitive.
3. Treatment time: 15 minutes.
4. Remarks: Since interferential current is an alternating current, the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: 3 min. MF = 10.000 Hz and AMF = 100 Hz.
12 min. MF = 10.000 Hz, AMF = 80 Hz, frequency sweep = 80 Hz and sweep cycle = 12/12 (in seconds).
6. Sensation: Mild, rhythmically occurring tingling sensation.
7. Intensity level: Just perceptible to clearly perceptible (sensoric stimulation).

IV. SUBACUTE PHASE (MF)

1. General: Programme for Medium Frequency current (interferential current, bipolar mode) using varying parameters in order to prevent accommodation.
2. Indication: Acute and subacute pain with symptoms of muscular hypertonia.
3. Treatment time: 15 minutes.
4. Remarks: Since interferential current is an alternating current, the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: 3 min. MF = 6.000 Hz and AMF = 60 Hz.
12 min. MF = 6.000 Hz, AMF = 30 Hz, frequency sweep = 50 Hz and sweep cycle = 6/6 (in seconds).
6. Sensation: Clearly perceptible, varying tingling sensation.
7. Intensity level: Strong perceptible, but not painful (just below motoric threshold).

V. CHRONIC PHASE (MF)

1. General: Programme for Medium Frequency current (interferential current, bipolar mode) using varying parameters in order to prevent accommodation.
2. Indication: Chronic pain disorders.
3. Treatment time: 20 minutes.
4. Remarks: Since interferential current is an alternating current, the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: 5 min. MF = 4.000 Hz and AMF = 20 Hz.
15 min. MF = 4.000 Hz, AMF = 1 Hz, frequency sweep = 29 Hz and sweep cycle = 1/1 (in seconds).
6. Sensation: Fast intermittent, tingling sensations.
7. Intensity level: Strongly perceptible up to motoric level.

VI. ACUTE PHASE (TENS)

1. General: Application of a comfortable form of TENS current for pain relief with frequency modulations and a slow sweep in order to prevent accommodation.
2. Indication: Acute pain disorders or syndromes.
3. Treatment time: 20 minutes (can be prolonged if necessary).
4. Remarks: Since this form of TENS current is biphasic (alternating), the risk for electrochemical damage (often called a chemical burn) is minimal.
However, negative electrode is more stimulating than positive electrode.
5. Programme: 5 min. Tens BF.ASYM. = 80 μ s, frequency 80 Hz.
15 min. Tens BF.ASYM. = 80 μ s, frequency 80 Hz, sweep 20 Hz and sweep cycle = 6/6 (in seconds).
6. Sensation: Slightly prickling, tingling.
7. Intensity level: Just perceptible, comfortable.

VII. SUBACUTE PHASE (TENS)

1. General: Application of a comfortable form of TENS current for pain relief with frequency modulations and a slow sweep in order to prevent accommodation.
2. Indication: For pain relief of disorders at the transition acute-chronic.
3. Treatment time: 25 minutes (can even be prolonged if necessary)
4. Remarks: Since this form of TENS current is biphasic (alternating), the risk for electrochemical damage (often called a chemical burn) is minimal.
However, negative electrode is more stimulating than positive electrode.
5. Programme: 5 min. BF.ASYM. = 200 μ s, frequency 50Hz
20 min. BF.ASYM. = 150 μ s, frequency 5 Hz, sweep 95 Hz and sweep cycle = 1/30 (in seconds).
6. Sensation: Rough tingling. Vibration.
7. Intensity level: From just perceptible up to patients' tolerance level.

VIII. CHRONIC PHASE (TENS)

1. General: Programme for pain relief at a chronic stage. First part of the programme slightly aggressive followed by a more milder phase.
2. Indication: Chronic pain disorders or syndromes.
3. Treatment time: 30 minutes.
4. Remarks: Since this form of TENS current is biphasic (alternating), the risk for electrochemical damage (often called a chemical burn) is minimal.
However, negative electrode is more stimulating than positive electrode.
5. Programme: 10 min. BF.ASYM. = 200 μ s, frequency 5Hz
20 min. BF.ASYM. = 100 μ s, frequency 85 Hz, burst = 2 Hz.
6. Sensation: Rhythmic vibrations. Muscular activity.
7. Intensity level: Strongly perceptible. From motoric level up to patients' tolerance level.

IX. DEEP PAIN: ACUTE PHASE (MF) (Endomed 582 ID only)

1. General: Programme for pain relief by means of two medium frequency alternating currents, causing interference (four-polar interferential).
2. Indication: Treatment of deeply felt pain (e.g. joint pain) at an acute stage.
3. Treatment time: 20 minutes.
4. Remarks: Since interferential current is an alternating current, the risk for electrochemical damage (often called a chemical burn) is minimal.
The electrode pairs from each circuit are placed in such a way, that the two currents cross.
5. Programme: 20 min. IF-4 pole, MF = 10.000 Hz, AMF = 100 Hz, frequency sweep = 50 Hz and sweep cycle = 12/12 (in seconds).
6. Sensation: Comfortable, rhythmically alternating stimulation from mild (100Hz) to strong (150 Hz).
7. Intensity level: Just perceptible (sensoric stimulation).

X. DEEP PAIN: SUBACUTE PHASE (TENS) (Endomed 582 ID only)

1. General: Programme for pain relief by means of burst TENS.
2. Indication: Treatment of deeply felt pain at the transition acute-chronic.
3. Treatment time: 30 minutes.
4. Remarks: Since this form of TENS current is biphasic (alternating), the risk for electrochemical damage (often called a chemical burn) is minimal.
However, negative electrode is more stimulating than positive electrode.
5. Programme: 30 min. TENS, BF.ASYM. = 150 μ s, frequency 80 Hz, burst = 2 Hz.
6. Sensation: Rhythmic vibrations, well tolerable. Muscular activity.
7. Intensity level: Up to motoric level. Strongly perceptible.

XI. DEEP PAIN: SUBACUTE PHASE (MF) (Endomed 582 ID only)

1. General: Programme for pain relief by means of two medium frequency alternating currents, causing interference (four-polar interferential).
2. Indication: Treatment of deeply felt pain at the transition acute-chronic.
3. Treatment time: 20 minutes.
4. Remarks: Since interferential current is an alternating current, the risk for electrochemical damage (often called a chemical burn) is minimal.
The electrode pairs from each circuit are placed in such a way, that the two currents cross.
5. Programme: 20 min. IF-4 pole, MF = 6.000 Hz, AMF = 80 Hz, frequency sweep = 40Hz and sweep cycle = 6/6 (in seconds).
6. Sensation: Comfortable, rhythmically alternating stimulation.
7. Intensity level: Just perceptible (sensoric stimulation).

XII. MYALGIA (Endomed 582 ID only)

1. General: Programme for muscular pain relief by means of two medium frequency Alternating currents, causing interference (four-polar interferential).
2. Indication: Myalgia (muscular pain).
3. Treatment time: 15 minutes.
4. Remarks: Since interferential current is an alternating current, the risk for electrochemical damage (often called a chemical burn) is minimal.
The electrode pairs from each circuit are placed in such a way, that the two currents cross.
5. Programme: 15 min. IF-4 pole + vector scan, MF = 6.000 Hz, AMF = 100 Hz, frequency sweep = 80 Hz and sweep cycle = 6/6 (in seconds).
6. Sensation: Comfortable, rhythmically alternating stimulation.
7. Intensity level: Just perceptible (sensoric stimulation).

I. SUPERFICIAL CIRCULATION IMPROVEMENT (DIA.)

1. General: Programme for improvement of the blood flow in superficial tissues.
2. Indication: Circulatory disorders.
3. Treatment time: 6 minutes.
4. Remarks: Diadynamic current is a Direct Current. At higher current densities DC may cause tissue damage. Intensity should not exceed 0.2 mA/cm² of effective treatment area. Preferably use thick pads/sponges (min. 1 cm).
5. Programme: 2 min. DF
2 min. CP
automatic polarity reverse
2 min. CP
6. Sensation: Slightly prickling, followed by alternating light to moderate vibrations.
7. Intensity level: Strongly perceptible.

II. CIRCULATION IMPROVEMENT (TENS)

1. General: Programme for improvement of the blood flow by means of activation of the 'muscle pump' function.
2. Indication: Circulatory disorders.
3. Treatment time: 10 minutes.
4. Remarks: Since this form of TENS current is biphasic (alternating), the risk for electrochemical damage (often called a chemical burn) is minimal. No difference between the negative electrode and positive electrode.
5. Programme: 10 min. BF.SYM. = 150 µs , frequency = 15 Hz, contraction rhythm = 1/1/1/3 (in seconds).
6. Sensation: Sequence of fast contractions (fast ramp-up and ramp-down) followed by brief intervals.
7. Intensity level: Clearly visible muscle contractions.

III. CIRCULATION IMPROVEMENT (MF)

1. General: Programme for improvement of the blood flow by means of activation of the 'muscle pump' function. Due to the application of a medium frequency current, presumably a more profound action is obtained.
2. Indication: Circulatory disorders.
3. Treatment time: 10 minutes.
4. Remarks: Since interferential current is an alternating current, the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: 10 min. MF = 4000 Hz, AMF = 50 Hz, sweep 50 Hz and sweep cycle = 1/1 (in seconds). Incorporated in a contraction rhythm of 1/1/1/3 (seconds).
6. Sensation: Sequence of fast contractions (fast ramp-up and ramp-down) followed by brief intervals.
7. Intensity level: Clearly visible muscle contractions.

IV. INTERMITTENT CLAUDICATION (TENS) (Endomed 582 ID only)

1. General: Programme for improvement of the blood flow in the lower extremities in case of ischemia.
2. Indication: Intermittent claudication.
3. Treatment time: 30 minutes.
4. Remarks: Since this form of TENS current is biphasic (alternating), the risk for electrochemical damage (often called a chemical burn) is minimal. No difference between the negative electrode and positive electrode.
5. Programme: 30 min. BF.SYM. = 200 µs , frequency = 4 Hz.
6. Sensation: Fast, rhythmically vibrating stimulation.
7. Intensity level: Clearly perceptible (sensoric level).

I. SINGULAR PROGRAMME (TENS)

1. General: Programme for muscle relaxation by reducing the muscular tonicity.
2. Indication: Muscular hypertonia.
3. Treatment time: 5 minutes.
4. Remarks: Since this form of TENS current is biphasic (alternating), the risk for electrochemical damage (often called a chemical burn) is minimal.
No difference between the negative electrode and positive electrode.
5. Programme: 5 min. BF.SYM. = 200 μ s, frequency = 10 Hz, sweep = 15 Hz and sweep cycle = 6/6 (in seconds).
6. Sensation: Strong vibration.
7. Intensity level: Until the required sensation is reached.

II. SEQUENTIAL PROGRAMME (MF)

1. General: Multiple programme for muscle relaxation.
2. Indication: Muscular hypertonia.
3. Treatment time: 6 minutes.
4. Remarks: Since interferential current is an alternating current, the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: 0-90 sec: MF = 2500 Hz, AMF = 10 Hz
90-186 sec: MF = 2500 Hz, AMF = 10 Hz, sweep = 15 Hz and sweep cycle = 12/12 (in seconds)
187-282 sec: MF = 2500 Hz, AMF = 10 Hz, sweep = 30 Hz and sweep cycle = 12/12 (in seconds)
283-360 sec: MF = 2500 Hz, AMF = 10 Hz
6. Sensation: Strong vibration.
7. Intensity level: Until the required sensation is reached.

I. MUSCLE STIMULATION (FARADISM)

1. General: Classic protocol to evoke a tetanic contraction by means of multiple electrical pulses.
2. Indication: First stage of therapeutic stimulation of striated muscles after reinnervation.
3. Treatment time: 15 minutes.
4. Remarks: It is recommended to evaluate the state of innervation by means of S-D curves. Faradic current is a Direct Current. At higher current densities DC may cause tissue damage. Intensity should not exceed 0.2 mA/cm² of effective treatment area. Preferably use thick pads/sponges (min. 1 cm).
5. Programme: Faradic current: pulse time = 1 ms, pause time = 19 ms. Integrated in a contraction rhythm of 1/1/1/6 (seconds).
6. Sensation: Moderate to fast contractions (fast ramp-up and ramp-down) followed by relatively long intervals.
7. Intensity level: Adjust intensity level until a tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

II. MUSCLE STIMULATION (TENS)

1. General: Programme to evoke a comfortable interrupted tetanic contraction.
2. Indication: First stage of rehabilitation of striated muscles.
3. Treatment time: 12 minutes
4. Remarks: Since this form of TENS current is biphasic (alternating), the risk for electrochemical damage (often called a chemical burn) is minimal. No difference between the negative electrode and positive electrode.
5. Programme: 12 min. BF.SYM. = 200 µs, frequency = 35 Hz. Incorporated in a contraction rhythm of 3/2/1/6 (seconds).
6. Sensation: Contraction pattern characterised by a slow ramp-up and a fast ramp-down followed by intervals.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

III. MUSCLE STIMULATION: PHASIC FIBRES (INITIAL STAGE)

1. General: Programme to evoke a comfortable interrupted tetanic contraction by means of TENS current.
2. Indication: First stage of rehabilitation of striated muscles with the objective to influence the twitch speed of fast-twitch (= phasic) fibres.
3. Treatment time: 10 minutes
4. Remarks: Since this form of TENS current is biphasic (alternating), the risk for electrochemical damage (often called a chemical burn) is minimal. No difference between the negative electrode and positive electrode.
5. Programme: 10 min. BF.SYM. = 100 µs, frequency = 50 Hz. Incorporated in a contraction rhythm of 1/4/1/12 (seconds).
6. Sensation: Comfortable muscle contractions in an on/off ratio of 1:3. Fast ramp-up and ramp-down followed by long intervals.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

IV. MUSCLE STIMULATION: PHASIC FIBRES (FINAL STAGE)

1. General: Programme to evoke a comfortable interrupted tetanic contraction by means of interferential current. This programme can be regarded as a continuation of the previous programme.
2. Indication: Advanced stage of rehabilitation of striated muscles with the objective to influence the twitch speed of fast-twitch (= phasic) fibres.
3. Treatment time: 15 minutes
4. Remarks: Since interferential current is an alternating current, the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: 15 min. MF = 2500 Hz, AMF = 80 Hz.
Incorporated in a contraction rhythm of 2/2/1/12 (seconds).
6. Sensation: Gradually increasing stimulation (slow ramp-up and fast ramp-down) followed by a long interval.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

V. MUSCLE STIMULATION: TONIC FIBRES (INITIAL STAGE)

1. General: Programme to evoke a comfortable interrupted tetanic contraction by means of interferential current.
2. Indication: First stage of rehabilitation of striated muscles with the objective to influence the twitch speed of slow-twitch (= tonic) fibres.
3. Treatment time: 15 minutes
4. Remarks: Since interferential current is an alternating current, the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: 15 min. MF = 2500 Hz, AMF = 20 Hz.
Incorporated in a contraction rhythm of 2/5/1/16 (seconds).
6. Sensation: Gradually increasing stimulation. Relatively long hold time followed by long interval.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

VI. MUSCLE STIMULATION: TONIC FIBRES (FINAL STAGE)

1. General: Programme to evoke a comfortable interrupted tetanic contraction by means of interferential current. This programme can be regarded as a continuation of the previous programme.
2. Indication: Advanced stage of rehabilitation of striated muscles with the objective to influence the twitch speed of slow-twitch (= tonic) fibres.
3. Treatment time: 15 minutes
4. Remarks: Since this form of TENS current is biphasic (alternating), the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: 15 min. BF.SYM. = 300 μ s, frequency = 20 Hz.
Incorporated in a contraction rhythm of 2/3/1/6 (seconds).
6. Sensation: Gradually increasing stimulation (slow ramp-up and fast ramp-down) followed by an interval.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

I. NEUTRAL

1. General: General programme for muscle strengthening by means of biphasic symmetrical TENS current.
2. Indication: Intensive muscle training with the objective to increase muscle strength without affecting the muscle fibre composition.
3. Treatment time: 32 minutes.
4. Remarks: Since this form of TENS current is biphasic (alternating), the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: 8 sequences of each 2 min. BF.SYM. = 200 μ s, frequency = 45 Hz. Incorporated in a contraction rhythm of 1/4/1/9 (seconds). An interval time of 2 minutes is given after each sequence.
6. Sensation: Sequence of gradually increasing and decreasing stimulation (fast ramp-up and fast ramp-down) followed by a long interval time.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

II. SLOW-TWITCH FIBRES

1. General: Specific programme for muscle strengthening by means of biphasic symmetrical TENS current.
2. Indication: Intensive muscle training with the objective to increase muscle strength of especially the slow-twitch (=tonic) fibres.
3. Treatment time: 40 minutes.
4. Remarks: Since this form of TENS current is biphasic (alternating), the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: 10 sequences of each 2 min. BF.SYM. = 300 μ s, frequency = 30 Hz. Incorporated in a contraction rhythm of 1/4/1/6 (seconds). An interval time of 2 minutes is given after each sequence.
6. Sensation: Sequence of gradually increasing and decreasing stimulation (fast ramp-up and fast ramp-down) followed by a long interval time.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

III. FAST-TWITCH FIBRES

1. General: Specific programme for muscle strengthening by means of biphasic symmetrical TENS current.
2. Indication: Intensive muscle training with the objective to increase muscle strength of especially the fast-twitch (=phasic) fibres.
3. Treatment time: 27 minutes.
4. Remarks: Since this form of TENS current is biphasic (alternating), the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: 6 sequences of each 1½ min. BF.SYM. = 100 μ s, frequency = 65 Hz. Incorporated in a contraction rhythm of 1/4/1/12 (seconds). An interval time of 3 minutes is given after each sequence.
6. Sensation: Sequence of gradually increasing and decreasing stimulation (fast ramp-up and fast ramp-down) followed by a long interval time.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

IV. TOTAL TRAINING

1. General: Combined programme for non-specific muscle strengthening by means of alternating current. The training itself is started with and terminated by a muscle relaxation programme (comparable to warming-up and cooling-down).
2. Indication: Intensive muscle training with the objective to increase muscle strength.
3. Treatment time: 35 minutes.
4. Remarks: Alternating currents: the risk for electrochemical damage is minimal.
5. Programme: 90 sec. MF = 2500 Hz, AMF = 10 Hz, 8 sequences of each 2 min. BF.SYM. = 200 μ s, frequency = 45 Hz.
Contraction rhythm = 1/4/1/9 (seconds).
Inter-sequence interval: 2 min.
90 sec. MF = 2500 Hz, AMF = 10 Hz.
6. Sensation: Sequence of gradually increasing and decreasing stimulation (fast ramp-up and fast ramp-down) followed by a long interval time. Strong vibrations at the start and at the end of the programme.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

I. SINGULAR PROGRAMME (TENS)

1. General: Simple programme for non-specific muscle strengthening.
2. Indication: In conditions where an increment of muscle force is required for non-injured and/or fully recovered muscles.
3. Treatment time: 33 minutes.
4. Remarks: Since this form of TENS current is biphasic (alternating), the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: 6 sequences of each 144 sec. BF.SYM. = 400 μ s, frequency = 35 Hz. Incorporated in a contraction rhythm of 2/3/1/12 (seconds). An interval time of 3 minutes is given after each sequence.
6. Sensation: Sequence of gradually increasing and decreasing stimulation (slow ramp-up and fast ramp-down) followed by a long interval time.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation.

II. SEQUENTIAL PROGRAMME

1. General: Complete programme for non-specific muscle strengthening, composed of a 'warming-up', 2 training phases, multiple intervals and 'cooling-down'.
2. Indication: In conditions, where an increment of muscle force is required for non-injured and/or fully recovered muscles.
3. Treatment time: 22 minutes.
4. Remarks: Alternating currents: the risk for electrochemical damage is minimal.
5. Programme: 90 sec. MF = 2500 Hz, AMF = 10 Hz, followed by an interval of 1 min.
4 min. BF.SYM = 400 μ s, frequency = 20 Hz. Incorporated in a contraction rhythm of 2/5/1/8 (=15 contractions).
2 minutes rest. 4 sequences, each 2 min. BF.SYM. = 400 μ s, frequency = 20 Hz, sweep = 45 Hz and sweep cycle = 12/12 (in sec.), followed by interval of 1 min. (total = 12 min.). 90 sec. MF = 2500 Hz, AMF = 10 Hz.
6. Sensation: Sequence of gradually increasing and decreasing stimulation, interrupted by intervals of different duration. Strong vibrations during first and final 90 seconds.
7. Intensity level: First 90 seconds: adjust intensity until a strong vibrating sensation is felt. Then adjust intensity level (during contraction phase) until a smooth tetanic contraction occurs. Discontinue treatment as soon as contraction strength drops.

III. SEQUENTIAL PROGRAMME (PHASIC FIBRES)

1. General: Complete programme for specific muscle strengthening, composed of a 'warming-up', 2 training phases, multiple intervals and 'cooling-down'.
2. Indication: In conditions, where an increment of the strength of non-injured and/or fully recovered phasic (fast-twitch) muscles is required.
3. Treatment time: 22 minutes.
4. Remarks: Alternating currents: the risk for electrochemical damage is minimal.
5. Programme: 90 sec. MF = 2500 Hz, AMF = 10 Hz, followed by an interval of 1 min.
4 min. BF.SYM = 400 μ s, frequency = 50 Hz. Incorporated in a contraction rhythm of 1/4/1/14 (= 12 contractions).
Next: 2 minutes rest.
4 sequences, each 2 min. BF.SYM. = 400 μ s, frequency = 50 Hz, sweep = 100 Hz, sweep cycle = 1/1 (in sec.) and contraction rhythm = 1/4/1/6 (sec.).
Followed by an interval of 1 min. (total = 12 minutes).
90 sec. MF = 2500 Hz, AMF = 10 Hz.
6. Sensation: Sequences of relatively short, not gradually increasing/decreasing stimulation, interrupted by intervals of different duration. Strong vibrations during first and final part.
7. Intensity level: First 90 seconds: adjust intensity until a strong vibrating sensation is felt. Then adjust intensity level (during contraction phase) until a smooth tetanic contraction occurs. Discontinue treatment as soon as contraction strength drops.

IV. SEQUENTIAL PROGRAMME (TONIC FIBRES)

1. General: Complete programme for specific muscle strengthening, composed of a 'warming-up', 2 training phases, multiple intervals and 'cooling-down'.
2. Indication: In conditions, where an increment of the strength of non-injured and/or fully recovered tonic (slow-twitch) muscles is required.
3. Treatment time: 22 minutes.
4. Remarks: Alternating currents: the risk for electrochemical damage is minimal.
5. Programme: 90 sec. MF = 2500 Hz, AMF = 10 Hz, followed by an interval of 1 min.
4 min. BF.SYM = 400 μ s, frequency = 20 Hz. Incorporated in a contraction rhythm of 1/6/1/8 (=16 contractions).
Next: 2 minutes rest.

4 sequences, each 2 min. BF.SYM. = 400 μ s, frequency = 10 Hz, sweep = 25 Hz, sweep cycle = 12/12 (in sec.), followed by interval of 1 min. (total = 12 min.).
90 sec. MF = 2500 Hz, AMF = 10 Hz.
6. Sensation: Sequences of relatively short, not gradually increasing/decreasing stimulation, followed by sequences with alternating rising and falling stimulation. Interrupted by intervals of different duration. Strong vibrations during first and final part.
7. Intensity level: First 90 seconds: adjust intensity until a strong vibrating sensation is felt. Then adjust intensity level (during contraction phase) until a smooth tetanic contraction occurs. Discontinue treatment as soon as contraction strength drops.

V. SINGULAR PROGRAMME (MF)

1. General: Simple programme for non-specific muscle strengthening.
2. Indication: In conditions where an increment of muscle force is required for non-injured and/or fully recovered muscles.
3. Treatment time: 15 minutes
4. Remarks: Since interferential current is an alternating current, the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: 10 sequences of each 1 min. MF = 2500 Hz, AMF = 20 Hz and sweep = 65 Hz. Incorporated in a contraction rhythm of 1/4/1/14 (seconds). Interval after each sequence: 30 seconds (total = 15 minutes).
6. Sensation: Sequence of gradually increasing and decreasing stimulation (fast ramp-up and fast ramp-down) followed by a long interval time.
7. Intensity level: Adjust intensity level (during contraction phase) until a smooth tetanic contraction occurs. Do not increase intensity any further during operation.

VI. PROGRAMME BASED ON 'RUSSIAN STIMULATION'

1. General: Intensive programme for muscle strengthening by means of interrupted medium frequency current. The training itself is started with and terminated by a muscle relaxation programme ('warming-up' and 'cooling-down').
2. Indication: In conditions where an increment of muscle force is required for healthy or fully recovered muscles.
3. Treatment time: 23 minutes.
4. Remarks: As 'Russian current' is a medium frequency alternating current, the risk for electrochemical damage is minimal.
5. Programme: 90 sec. MF = 2500 Hz, frequency = 10 Hz. 8 sequences of each 1½ min. MF = 2500 Hz, frequency = 50 Hz, and ratio 1:1.
Incorporated in a contraction rhythm of 1/7/1/9 (seconds). Each sequence is followed by an interval of 1 minute.
90 sec. MF = 2500 Hz, AMF = 10 Hz.
6. Sensation: Sequences of fast, powerful contractions, interrupted by intervals. Strong vibrations at the start and at the end of the programme.
7. Intensity level: First 90 seconds: adjust intensity until a strong vibrating sensation is felt. Then adjust intensity level (during contraction phase) until patient's tolerance level. Discontinue treatment if contraction strength weakens.

I. BURSITIS (MF)

1. General: Programme for the improvement of the blood flow by means of medium frequency current (bipolar interferential), as part of a treatment protocol to reduce the inflammatory state.
2. Indication: Subacute and chronic bursitis, e.g. bursitis prae-patellaris.
3. Treatment time: 14 minutes.
4. Remarks: Since interferential current is an alternating current, the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: 14 min. bipolar interferential: MF = 2500 Hz, AMF = 80 Hz, sweep 50 Hz and sweep cycle = 6/6 (in seconds).
6. Sensation: Strong, continuous stimulation.
7. Intensity level: Clearly perceptible up to strongly perceptible.

II. CAPSULITIS (MF)

1. General: Programme for the improvement of the blood flow by means of medium frequency currents (four polar interferential), as part of a treatment protocol to reduce the inflammatory state of the articular capsule.
2. Indication: Subacute and chronic capsulitis, e.g. capsulitis of the shoulder joint.
3. Treatment time: 15 minutes.
4. Remarks: Since interferential current is an alternating current, the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: 15 min. four polar interferential: MF = 4000 Hz, AMF = 100 Hz, sweep 80 Hz and sweep cycle = 6/6 (in seconds).
6. Sensation: Tolerable, fast rhythmic stimulation.
7. Intensity level: Clearly perceptible up to strongly perceptible.

III. EPICONDYLITIS (TENS)

1. General: Programme for the improvement of the blood flow by means of TENS current, as part of a treatment protocol to reduce the symptoms of a 'tennis elbow'.
2. Indication: Subacute and chronic tendinitis and/or periostitis of the epicondylus, e.g. epicondylitis lateralis humeri.
3. Treatment time: 15 minutes.
4. Remarks: Since this form of TENS current is biphasic (alternating), the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: TENS BF.SYM. = 75 μ s, frequency = 80 Hz, sweep 20 Hz and sweep cycle = 6/6 (in seconds).
6. Sensation: Comfortable, fast rhythmic stimulation.
7. Intensity level: Strongly perceptible up to patient's tolerance level.

IV. TENDINITIS (MF)

1. General: Programme for the improvement of the blood flow by means of a medium frequency current (bipolar interferential), as part of a treatment protocol to reduce the inflammatory state of a tendon.
2. Indication: Subacute and chronic tendinitis.
3. Treatment time: 14 minutes.
4. Remarks: Since interferential current is an alternating current, the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: 14 min. bipolar interferential: MF = 2500 Hz, AMF = 80 Hz, sweep 80 Hz and sweep cycle = 1/30/1/30 (in seconds).
6. Sensation: Intruding, strong vibrations.
7. Intensity level: Clearly perceptible up to strongly perceptible.

I. STRESS INCONTINENCE

1. General: Programme for improving the strength of the striated pelvic floor muscles by means of a comfortable TENS current.
2. Indication: Stress incontinence.
3. Treatment time: 15 minutes.
4. Remarks: Since this form of TENS current is biphasic (alternating), the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: TENS BF.SYM. = 150 μ s, frequency = 35 Hz, incorporated in a contraction rhythm of 7/25 (seconds) during 15 minutes.
6. Sensation: Comfortable, but clearly perceptible contraction of the pelvic floor, followed by relaxation.
7. Intensity level: Adjust intensity level (during contraction phase) until a smooth, clearly perceptible contraction of the pelvic floor occurs. The intensity is not further increased during operation.

II. URGE INCONTINENCE

1. General: Programme for the inhibition of the bladder contraction and/or normalization of the tone by means of a comfortable TENS current.
2. Indication: Urge incontinence.
3. Treatment time: 15 minutes.
4. Remarks: Since this form of TENS current is biphasic (alternating), the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: TENS BF.SYM. = 100 μ s, frequency = 5 Hz, continuously during 15 min.
6. Sensation: Comfortable tingling in a low beat frequency.
7. Intensity level: Adjust the intensity level until clearly perceptible (sensoric stimulation).

III. MIXED INCONTINENCE

1. General: Programme, in which successively pelvic floor muscles are activated and the bladder contraction is inhibited.
2. Indication: Mixed incontinence.
3. Treatment time: 25 minutes.
4. Remarks: Since this form of TENS current is biphasic (alternating), the risk for electrochemical damage (often called a chemical burn) is minimal.
5. Programme: TENS BF.SYM. = 150 μ s, frequency = 35 Hz, incorporated in a contraction rhythm of 7/25 (seconds) during 15 minutes followed by: TENS BF.SYM. = 100 μ s, frequency = 5 Hz, continuously during 10 min.
6. Sensation: Clearly perceptible contraction followed by a comfortable tingling sensation.
7. Intensity level: First adjust the intensity level a clearly perceptible contraction of the pelvic floor occurs (motoric level). The intensity during the second phase should be set at clearly perceptible (sensoric stimulation).

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