
# Sonopuls 492 Therapy guide Protocols

Protocollen  
Protocols  
Protokolle  
Protocole  
Protocolos

Copyright:

**Enraf  
Nonius**

Enraf-Nonius B.V.  
P.O. Box 12080  
3004 GB ROTTERDAM  
The Netherlands  
Tel.: +31(0)10 - 20 30 600  
Fax: +31(0)10 - 20 30 699

Article number : 1497.761-42  
December 2005

# **Sonopuls 492**

## **Therapy guide**

### **Protocols**

**Protocollen**  
**Protocols**  
**Protokolle**  
**Protocole**  
**Protocolos**

Nederlands .....	5
English .....	25
Deutsch .....	45
français .....	65
Español .....	85

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The Sonopuls 492 offers the opportunity to store frequently used settings in a memory. Single programs can be stored under program numbers 0 till 9. Sequential programs can be stored under cluster number 1. The pre-programmed programs can be retrieved as from cluster number 2.

<b>2. PAIN RELIEF .....</b>	<b>27</b>
2.1. SUPERFICIAL PAIN (DIA.) .....	27
2.2. NEUROGENIC (DIA.) .....	27
2.3. ACUTE PHASE (MF) .....	27
2.4. SUBACUTE PHASE (MF) .....	27
2.5. CHRONIC PHASE (MF) .....	28
2.6. ACUTE PHASE (TENS) .....	28
2.7. SUBACUTE PHASE (TENS) .....	28
2.8. CHRONIC PHASE (TENS) .....	28
2.9. DEEP PAIN: SUBACUTE PHASE (TENS) .....	28
2.10. POSTHERPETIC NEURALGIA DEX. TH5-TH7 (TENS) .....	29
2.11. DEEP PAIN: ACUTE PHASE (MF) .....	29
2.12. DEEP PAIN: SUBACUTE PHASE (MF) .....	29
2.13. MYALGIA (MF) .....	29
<b>3. CIRCULATION IMPROVEMENT .....</b>	<b>30</b>
3.1. SUPERFICIAL CIRCULATION IMPROVEMENT (DIA.) .....	30
3.2. CIRCULATION IMPROVEMENT (TENS) .....	30
3.3. CIRCULATION IMPROVEMENT (MF) .....	30
3.4. INTERMITTENT CLAUDICATION (TENS) .....	30
<b>4. MUSCLE RELAXATION - REDUCTION OF TONICITY .....</b>	<b>31</b>
4.1. SINGULAR PROGRAMME (TENS) .....	31
4.2. SEQUENTIAL PROGRAMME (MF) .....	31
<b>5. NERVE LESIONS .....</b>	<b>32</b>
5.1. MUSCLE STIMULATION (FARADISM) .....	32
5.2. MUSCLE STIMULATION (TENS) .....	32
5.3. MUSCLE STIMULATION: PHASIC FIBRES (INITIAL STAGE) (TENS) .....	32
5.4. MUSCLE STIMULATION: PHASIC FIBRES (FINAL STAGE) (MF) .....	32
5.5. MUSCLE STIMULATION: TONIC FIBRES (INITIAL STAGE) (MF) .....	33
5.6. MUSCLE STIMULATION: TONIC FIBRES (FINAL STAGE) (TENS) .....	33
<b>6. MUSCLE TRAINING .....</b>	<b>34</b>
6.1. NEUTRAL (TENS) .....	34
6.2. SLOW-TWITCH FIBRES (TENS) .....	34
6.3. FAST-TWITCH FIBRES (TENS) .....	34
6.4. TOTAL TRAINING (MF,TENS) .....	34
6.5. POST-OPERATIVE (MF) .....	35
6.6. ENDURANCE (AEROBIC) (TENS) .....	35
6.7. ENDURANCE (MAX. AEROBIC) (TENS) .....	35
6.8. MAXIMAL CONTRACTION (TENS) .....	35
<b>7. MUSCLE STRENGTHENING .....</b>	<b>36</b>
7.1. SINGULAR PROGRAMME (TENS) .....	36
7.2. SEQUENTIAL MUSCLE STRENGTHENING PROGRAMME (MF,TENS) .....	36
7.3. SEQUENTIAL PROGRAMME (PHASIC FIBRES) (MF,TENS) .....	36
7.4. SEQUENTIAL PROGRAMME (TONIC FIBRES) (MF,TENS) .....	37
7.5. SINGULAR PROGRAMME (MF) .....	37
7.6. PROGRAMME BASED ON 'RUSSIAN STIMULATION' .....	37

<b>8. INFLAMMATION .....</b>	<b>38</b>
8.1. BURSITIS (MF) .....	38
8.2. EPICONDYLITIS (TENS) .....	38
8.3. TENDINITIS (MF) .....	38
8.4. CAPSULITIS (MF) .....	38
<b>9. URINARY INCONTINENCE .....</b>	<b>39</b>
9.1. STRESS INCONTINENCE (TENS) .....	39
9.2. URGE INCONTINENCE (TENS) .....	39
9.3. MIXED INCONTINENCE (TENS) .....	39
<b>10. TRANSVERSAL LESION OF THE SPINAL CORD AND C.V.A. ....</b>	<b>40</b>
10.1. TROPHIC IMPROVEMENT (TENS) .....	40
10.2. CONTRACTURE RESTRICTION (TENS) .....	40
10.3. TONE NORMALIZATION HYPERTONIA (TENS) .....	40
<b>11. AMYOTROFIE .....</b>	<b>41</b>
11.1. UPPER EXTREMITY AND TRUNK (TENS) .....	41
11.2. LOWER EXTREMITY (TENS) .....	41
<b>13. ULTRASOUND .....</b>	<b>42</b>
13.1. ARTHRITIS .....	42
13.2. BURSITIS/CAPSULITIS .....	42
13.3. EPICONDYLITIS .....	42
13.4. CONTUSIONS/DISTORSIONS .....	42
13.5. LIGAMENT SPRAIN .....	42
13.6. SCAR TISSUE .....	42
13.7. MUSCLE STRAIN .....	42
13.8. SPONDYLOSIS DEFORMANS .....	42
13.9. TENDINITIS .....	43

**2.1. SUPERFICIAL PAIN (DIA.)**

1. General: Programme for the treatment of pain related disorders by using the milder current forms of diadynamic current as described by dr. Bernard.
2. Indication: Superficial, sharp and/or dull pain.
3. Treatment time: 8 minutes.
4. Remarks: Diadynamic current is a Direct Current. At higher current densities DC may cause tissue damage. Intensity should not exceed 0.2mA/cm<sup>2</sup> of effective treatment area. Preferably use thick pads/sponges (min. 1 cm).
5. Programme: 2 min. DF, 3 min. LP, automatic polarity reverse 3 min. LP
6. Sensation: Tickling, prickling followed by alternately light to moderate vibrations.
7. Intensity level: Clearly to strongly perceptible, preferably just below motoric threshold.

**2.2. NEUROGENIC (DIA.)**

1. General: Programme for the treatment of neurogenic pain by using the CP form of diadynamic current (CP = courant modulé en Courtes Périodes) and automatic polarity reverse.
2. Indication: Neuralgia, causalgia, post-herpetic neuralgia, phantom limb pain (besides pain relief, local circulation is increased).
3. Treatment time: 6 minutes.
4. Remarks: Diadynamic current is a Direct Current. At higher current densities DC may cause tissue damage. Intensity should not exceed 0.2 mA/cm<sup>2</sup> of effective treatment area. Preferably use thick pads/sponges (min. 1cm).
5. Programme: 3 min. CP, automatic polarity reverse, 3 min. CP
6. Sensation: Alternately light and rough vibrations.
7. Intensity level: Clearly perceptible, not painful. Preferably just below motoric threshold.

**2.3. ACUTE PHASE (MF)**

1. General: Programme for Medium Frequency current (interferential current, bipolar mode) using varying parameters in order to prevent accommodation.
2. Indication: Acute pain in absence of muscular hypertonia or as initial treatment for patients who are extremely sensitive.
3. Treatment time: 15 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 10 kHz, AMF 100 Hz, 3 min.  
IF, 10 kHz, AMF 80 Hz, frequency sweep 80 Hz, sweep cycle 12/12, 15 min.
6. Sensation: Mild, rhythmically occurring tingling sensation.
7. Intensity level: Just perceptible to clearly perceptible (sensoric stimulation).

**2.4. SUBACUTE PHASE (MF)**

1. General: Programme for Medium Frequency current (interferential current, bipolar mode) using varying parameters in order to prevent accommodation.
2. Indication: Acute and subacute pain with symptoms of muscular hypertonia.
3. Treatment time: 15 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 6 kHz, AMF 60 Hz, 3 min.  
IF, 6 kHz, AMF 30 Hz, frequency sweep 50 Hz, sweep cycle 6/6, 12 min.
6. Sensation: Clearly perceptible, varying tingling sensation.
7. Intensity level: Strong perceptible, but not painful (just below motoric threshold).

## 2.5. CHRONIC PHASE (MF)

1. General: Programme for Medium Frequency current (interferential current, bipolar mode) using varying parameters in order to prevent accommodation.
2. Indication: Chronic pain disorders.
3. Treatment time: 20 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 4 kHz, AMF 20 Hz, 5 min.  
IF, 4 kHz, AMF 1 Hz, frequency sweep 29 Hz, sweep cycle 1/1, 15 min.
6. Sensation: Fast intermittent, tingling sensations.
7. Intensity level: Strongly perceptible up to motoric level.

## 2.6. ACUTE PHASE (TENS)

1. General: Application of a comfortable form of TENS current for pain relief with frequency modulations and a slow sweep in order to prevent accommodation.
2. Indication: Acute pain disorders or syndromes.
3. Treatment time: 20 minutes (can be prolonged if necessary).
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.ASYM., 75  $\mu$ s, 80 Hz, 5 min.  
BF.ASYM., 75  $\mu$ s, 80 Hz, frequency sweep 20 Hz, sweep cycle 6/6, 15 min.
6. Sensation: Slightly prickling, tingling.
7. Intensity level: Just perceptible, comfortable.

## 2.7. SUBACUTE PHASE (TENS)

1. General: Application of a comfortable form of TENS current for pain relief with frequency modulations and a slow sweep in order to prevent accommodation.
2. Indication: For pain relief of disorders at the transition acute-chronic.
3. Treatment time: 25 minutes (can even be prolonged if necessary)
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.ASYM., 200  $\mu$ s, 50 Hz, 5 min.  
BF.ASYM., 150  $\mu$ s, 5 Hz, frequency sweep 90 Hz, sweep cycle 1/30, 20 min.
6. Sensation: Rough tingling vibration.
7. Intensity level: From just perceptible up to patients tolerance level.

## 2.8. CHRONIC PHASE (TENS)

1. General: Programme for pain relief at a chronic stage. First part of the programme slightly aggressive followed by a more milder phase.
2. Indication: Chronic pain disorders or syndromes.
3. Treatment time: 30 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.ASYM., 200  $\mu$ s, 5 Hz, 10 min.  
BF.ASYM., 100  $\mu$ s, 80 Hz, burst 2 Hz, 20 min.
6. Sensation: Rhythmic vibrations. Muscular activity.
7. Intensity level: Strongly perceptible. From motoric level up to patients tolerance level.

## 2.9. DEEP PAIN: SUBACUTE PHASE (TENS)

1. General: Programme for pain relief by means of burst TENS.
2. Indication: Treatment of deeply felt pain at the transition acute-chronic.
3. Treatment time: 30 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.ASYM., 150  $\mu$ s, 80 Hz, burst 2 Hz, 30 min.
6. Sensation: Rhythmic vibrations, well tolerable. Muscular activity.
7. Intensity level: Up to motoric level. Strongly perceptible.

**2.10. POSTHERPETIC NEURALGIA DEX. TH5-TH7 (TENS)**

1. General: Application of a comfortable form of TENS current for pain relief.
2. Indication: Pain in the affected segments, which is still present 30 day after eruption of the vesicles.
3. Treatment time: 30 min.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.ASYM., 75  $\mu$ s, 100 Hz, 30 min.
6. Sensation: Prickling tingling.
7. Intensity level: Up to patients tolerance level.

**2.11. DEEP PAIN: ACUTE PHASE (MF)**

1. General: Programme for pain relief by means of two medium frequency alternating currents, causing interference (four-polar interferential).
2. Indication: Treatment of deeply felt pain (e.g. joint pain) at an acute stage.
3. Treatment time: 20 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 4p, 10 kHz, AMF 100 Hz, frequency sweep 50 Hz, sweep cycle 6/6, 20 min.
6. Sensation: Comfortable, rhythmically alternating stimulation from mild (100Hz) to strong (150 Hz).
7. Intensity level: Just perceptible (sensoric stimulation).

**2.12. DEEP PAIN: SUBACUTE PHASE (MF)**

1. General: Programme for pain relief by means of two medium frequency alternating currents, causing interference (four-polar interferential).
2. Indication: Treatment of deeply felt pain at the transition acute-chronic.
3. Treatment time: 20 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 4p, 6 kHz, AMF 80 Hz, frequency sweep 40 Hz, sweep cycle 6/6, 20 min.
6. Sensation: Comfortable, rhythmically alternating stimulation.
7. Intensity level: Just perceptible (sensoric stimulation).

**2.13. MYALGIA (MF)**

1. General: Programme for muscular pain relief by means of two medium frequency Alternating currents, causing interference (four-polar interferential).
2. Indication: Myalgia (muscular pain).
3. Treatment time: 15 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 4p, isoplanary vector field, 6 kHz, AMF 100 Hz, frequency sweep 80 Hz, sweep cycle 6/6, 15 min.
6. Sensation: Comfortable, rhythmically alternating stimulation.
7. Intensity level: Just perceptible (sensoric stimulation).

**3.1. SUPERFICIAL CIRCULATION IMPROVEMENT (DIA.)**

1. General: Programme for improvement of the blood flow in superficial tissues.
2. Indication: Circulatory disorders.
3. Treatment time: 6 minutes.
4. Remarks: Diadynamic current is a Direct Current. At higher current densities DC may cause tissue damage. Intensity should not exceed 0.2 mA/cm<sup>2</sup> of effective treatment area. Preferably use thick pads/sponges (min. 1 cm).
5. Programme: 2 min. DF, 2 min. CP, automatic polarity reverse, 2 min. CP
6. Sensation: Slightly prickling, followed by alternating light to moderate vibrations.
7. Intensity level: Strongly perceptible.

**3.2. CIRCULATION IMPROVEMENT (TENS)**

1. General: Programme for improvement of the blood flow by means of activation of the muscle pump function.
2. Indication: Circulatory disorders.
3. Treatment time: 10 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 150  $\mu$ s, 15 Hz, contraction rhythm 1/1/1/3, 10 min.
6. Sensation: Sequence of fast contractions (fast ramp-up and ramp-down) followed by brief intervals.
7. Intensity level: Clearly visible muscle contractions.

**3.3. CIRCULATION IMPROVEMENT (MF)**

1. General: Programme for improvement of the blood flow by means of activation of the muscle pump function. Due to the application of a medium frequency current, presumably a more profound action is obtained.
2. Indication: Circulatory disorders.
3. Treatment time: 10 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 4 kHz, AMF 50 Hz, frequency sweep 50 Hz, sweep cycle 1/1, contraction rhythm 1/1/1/3, 10 min.
6. Sensatio: Sequence of fast contractions (fast ramp-up and ramp-down) followed by brief intervals.
7. Intensity level: Clearly visible muscle contractions.

**3.4. INTERMITTENT CLAUDICATION (TENS)**

1. General: Programme for improvement of the blood flow in the lower extremities in case of ischemia.
2. Indication: Intermittent claudication.
3. Treatment time: 30 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 200  $\mu$ s, 4 Hz, 30 min.
6. Sensation: Fast, rhythmically vibrating stimulation.
7. Intensity level: Clearly perceptible (sensoric level).

**4.1. SINGULAR PROGRAMME (TENS)**

1. General: Programme for muscle relaxation by reducing the muscular tonicity.
2. Indication: Muscular hypertonia.
3. Treatment time: 5 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 200  $\mu$ s, 10 Hz, frequency sweep 15 Hz, sweep cycle 6/6, 5 min.
6. Sensation: Strong vibration.
7. Intensity level: Until the required sensation is reached.

**4.2. SEQUENTIAL PROGRAMME (MF)**

1. General: Multiple programme for muscle relaxation.
2. Indication: Muscular hypertonia.
3. Treatment time: 6 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme:  
IF, 2,5 kHz, AMF 10 Hz, 1,5 min.  
IF, 2,5 kHz, AMF 10 Hz, frequency sweep 15 Hz, sweep cycle 12/12, 1,5 min.  
IF, 2,5 kHz, AMF 10 Hz, frequency sweep 30 Hz, sweep cycle 12/12, 1,5 min.  
IF, 2,5 kHz, AMF 10 Hz, 1,5 min.
6. Sensation: Strong vibration.
7. Intensity level: Until the required sensation is reached.

**5.1. MUSCLE STIMULATION (FARADISM)**

1. General: Classic protocol to evoke a tetanic contraction by means of multiple electrical pulses.
2. Indication: First stage of therapeutic stimulation of striated muscles after reinnervation.
3. Treatment time: 15 minutes.
4. Remarks: It is recommended to evaluate the state of innervation by means of S-D curves. Faradic current is a Direct Current. At higher current densities DC may cause tissue damage. Intensity should not exceed 0.2 mA/cm<sup>2</sup> of effective treatment area. Preferably use thick pads/sponges (min. 1 cm).
5. Programme: Faradic current: pulse time 1 ms, pause time 20 ms, contraction rhythm of 1/1/1/6, 15 min.
6. Sensation: Moderate to fast contractions (fast ramp-up and ramp-down) followed by relatively long intervals.
7. Intensity level: Adjust intensity level until a tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

**5.2. MUSCLE STIMULATION (TENS)**

1. General: Programme to evoke a comfortable interrupted tetanic contraction.
2. Indication: First stage of rehabilitation of striated muscles.
3. Treatment time: 12 minutes
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 200  $\mu$ s, 35 Hz, contraction rhythm 3/2/1/6, 12 min.
6. Sensation: Contraction pattern characterised by a slow ramp-up and a fast ramp-down followed by intervals.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

**5.3. MUSCLE STIMULATION: PHASIC FIBRES (INITIAL STAGE) (TENS)**

1. General: Programme to evoke a comfortable interrupted tetanic contraction by means of TENS current.
2. Indication: First stage of rehabilitation of striated muscles with the objective to influence the twitch speed of fast-twitch (= phasic) fibres.
3. Treatment time: 10 minutes
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 100  $\mu$ s, 50 Hz, contraction rhythm 1/4/1/12, 10 min.
6. Sensation: Comfortable muscle contractions in an on/off ratio of 1:3. Fast ramp-up and ramp-down followed by long intervals.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

**5.4. MUSCLE STIMULATION: PHASIC FIBRES (FINAL STAGE) (MF)**

1. General: Programme to evoke a comfortable interrupted tetanic contraction by means of interferential current. This programme can be regarded as a continuation of programme 5.3.
2. Indication: Advanced stage of rehabilitation of striated muscles with the objective to influence the twitch speed of fast-twitch (= phasic) fibres.
3. Treatment time: 15 minutes
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 2,5 kHz, AMF 80 Hz, contraction rhythm 2/2/1/12, 15 min.
6. Sensation: Gradually increasing stimulation (slow ramp-up and fast ramp-down) followed by a long interval.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

**5.5. MUSCLE STIMULATION: TONIC FIBRES (INITIAL STAGE) (MF)**

1. General: Programme to evoke a comfortable interrupted tetanic contraction by means of interferential current.
2. Indication: First stage of rehabilitation of striated muscles with the objective to influence the twitch speed of slow-twitch (= tonic) fibres.
3. Treatment time: 15 minutes
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 2,5 kHz, AMF 20 Hz, contraction rhythm 2/5/1/16, 15 min.
6. Sensation: Gradually increasing stimulation. Relatively long hold time followed by long interval.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

**5.6. MUSCLE STIMULATION: TONIC FIBRES (FINAL STAGE) (TENS)**

1. General: Programme to evoke a comfortable interrupted tetanic contraction by means of interferential current. This programme can be regarded as a continuation of programme 5.5.
2. Indication: Advanced stage of rehabilitation of striated muscles with the objective to influence the twitch speed of slow-twitch (= tonic) fibres.
3. Treatment time: 15 minutes
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 300  $\mu$ s, 20 Hz, contraction rhythm 2/3/1/6, 15 min.
6. Sensation: Gradually increasing stimulation (slow ramp-up and fast ramp-down) followed by an interval.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

### 6.1. NEUTRAL (TENS)

1. General: General programme for muscle strengthening by means of biphasic symmetrical TENS current.
2. Indication: Intensive muscle training with the objective to increase muscle strength without affecting the muscle fibre composition.
3. Treatment time: 32 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: 8 sequences of each 2 min. BF.SYM., 200  $\mu$ s, 45 Hz, contraction rhythm 1/4/1/9, interval time of 2 minutes after each sequence.
6. Sensation: Sequence of gradually increasing and decreasing stimulation (fast ramp-up and fast ramp-down) followed by a long interval time.
7. Intensity level: Adjust intensity level until the required contraction is reached. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

### 6.2. SLOW-TWITCH FIBRES (TENS)

1. General: Specific programme for muscle strengthening by means of biphasic symmetrical TENS current.
2. Indication: Intensive muscle training with the objective to increase muscle strength of especially the slow-twitch (=tonic) fibres.
3. Treatment time: 40 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: 10 sequences of each 2 min. BF.SYM., 300  $\mu$ s, 30 Hz, contraction rhythm 1/4/1/6, interval time of 2 minutes after each sequence.
6. Sensation: Sequence of gradually increasing and decreasing stimulation (fast ramp-up and fast ramp-down) followed by a long interval time.
7. Intensity level: Adjust intensity level until the required contraction is reached. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

### 6.3. FAST-TWITCH FIBRES (TENS)

1. General: Specific programme for muscle strengthening by means of biphasic symmetrical TENS current.
2. Indication: Intensive muscle training with the objective to increase muscle strength of especially the fast-twitch (=phasic) fibres.
3. Treatment time: 27 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: 6 sequences of each 1½ min. BF.SYM., 100  $\mu$ s, 60 Hz, contraction rhythm 1/4/1/12, interval time of 3 minutes after each sequence.
6. Sensation: Sequence of gradually increasing and decreasing stimulation (fast ramp-up and fast ramp-down) followed by a long interval time.
7. Intensity level: Adjust intensity level until the required contraction is reached. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

### 6.4. TOTAL TRAINING (MF,TENS)

1. General: Combined programme for non-specific muscle strengthening by means of alternating current. The training itself is started with and terminated by a muscle relaxation programme (comparable to warming-up and cooling-down).
2. Indication: Intensive muscle training with the objective to increase muscle strength.
3. Treatment time: 35 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 2,5 kHz, AMF 10 Hz, 1,5 min.  
8 sequences of each 2 min. BF.SYM., 200  $\mu$ s, 45 Hz, contraction rhythm 1/4/1/9, interval time of 2 minutes after each sequence.  
IF, 2,5 kHz, AMF 10 Hz, 1,5 min.
6. Sensation: Sequence of gradually increasing and decreasing stimulation (fast ramp-up and fast ramp-down) followed by a long interval time. Strong vibrations at the start and at the end of the programme.
7. Intensity level: Adjust intensity level until the required contraction is reached. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

### 6.5. POST-OPERATIVE (MF)

1. General: General programme for muscle strengthening post-operative at which no tetanic contraction occurs.
2. Indication: First stage of muscle rehabilitation with the objective to increase muscle strength without affecting the muscle fibre composition.
3. Treatment time: 7 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 2 kHz, AMF 20 Hz, 7 min.
6. Sensation: Perceptible contraction.
7. Intensity level: Adjust intensity level until the required contraction is reached. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

### 6.6. ENDURANCE (AEROBIC) (TENS)

1. General: Specific programme for aerobic endurance training by means of biphasic symmetrical TENS current.
2. Indication: Muscle rehabilitation with the objective to increase the endurance of the muscle.
3. Treatment time: 20 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 200  $\mu$ s, 35 Hz, contraction rhythm 1/5/0/6, 20 min.
6. Sensation: Gradually increasing stimulation followed by an interval.
7. Intensity level: Adjust intensity level until the required contraction is reached. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

### 6.7. ENDURANCE (MAX. AEROBIC) (TENS)

1. General: Specific programme for maximal aerobic endurance training by means of biphasic symmetrical TENS current.
2. Indication: Muscle rehabilitation with the objective to increase the maximal endurance of the muscle.
3. Treatment time: 10 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 200  $\mu$ s, 50 Hz, contraction rhythm 1/5/0/15, 10 min.
6. Sensation: Gradually increasing stimulation followed by an interval.
7. Intensity level: Adjust intensity level until the required contraction is reached. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

### 6.8. MAXIMAL CONTRACTION (TENS)

1. General: Specific programme for training of the maximal contraction of a muscle (-group).
2. Indication: Intensive muscle training with the objective to increase the maximal contraction of a muscle.
3. Treatment time: 7 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 200  $\mu$ s, 70 Hz, contraction rhythm 1/5/0/25, 7 min.
6. Sensation: Gradually increasing stimulation followed by an interval.
7. Intensity level: Adjust intensity level until the required contraction is reached. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

### 7.1. SINGULAR PROGRAMME (TENS)

1. General: Simple programme for non-specific muscle strengthening.
2. Indication: In conditions where an increment of muscle force is required for non-injured and/or fully recovered muscles.
3. Treatment time: 33 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: 6 sequences of each 150 sec. BF.SYM., 400  $\mu$ s, 35 Hz, contraction rhythm 2/3/1/12, interval time of 3 minutes after each sequence.
6. Sensation: Sequence of gradually increasing and decreasing stimulation (slow ramp-up and fast ramp-down) followed by a long interval time.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation.

### 7.2. SEQUENTIAL MUSCLE STRENGTHENING PROGRAMME (MF,TENS)

1. General: Complete programme for non-specific muscle strengthening, composed of a 'warming-up', 2 training phases, multiple intervals and 'cooling-down'.
2. Indication: In conditions, where an increment of muscle force is required for non-injured and/or fully recovered muscles.
3. Treatment time: 22 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 2,5 kHz, AMF 10 Hz, followed by an interval of 1 min.  
BF.SYM., 400  $\mu$ s, 20 Hz, contraction rhythm 2/5/1/8, 4 min, followed by an interval of 2 min.  
4 sequences of each 2 min. BF.SYM., 400  $\mu$ s, 20 Hz, frequency sweep 45 Hz, sweep cycle 12/12, followed by an interval of 1 min.  
IF, 2,5 kHz, AMF 10 Hz 1,5 min.
6. Sensation: Sequence of gradually increasing and decreasing stimulation, interrupted by intervals of different duration. Strong vibrations during first and final 90 seconds.
7. Intensity level: First 90 seconds: adjust intensity until a strong vibrating sensation is felt. Than adjust intensity level (during contraction phase) until a smooth tetanic contraction occurs. Discontinue treatment as soon as contraction strength drops.

### 7.3. SEQUENTIAL PROGRAMME (PHASIC FIBRES) (MF,TENS)

1. General: Complete programme for specific muscle strengthening, composed of a 'warming-up', 2 training phases, multiple intervals and 'cooling-down'.
2. Indication: In conditions, where an increment of the strength of non-injured and/or fully recovered phasic (fast-twitch) muscles is required.
3. Treatment time: 22 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 2p, 2,5 kHz, AMF 10 Hz, followed by an interval of 1 min.  
BF.SYM., 400  $\mu$ s, 50 Hz, contraction rhythm 1/4/1/14, 4 min. followed by an interval of 2 min.  
4 sequences of each 2 min. BF.SYM., 400  $\mu$ s, 50 Hz, frequency sweep 100 Hz, sweep cycle 1/1, contraction rhythm 1/4/1/6, followed by an interval of 1 min.  
IF, 2p, 2,5 kHz, AMF 10 Hz, 1,5 min.
6. Sensation: Sequences of relatively short, not gradually increasing/decreasing stimulation, interrupted by intervals of different duration. Strong vibrations during first and final part.
7. Intensity level: First 90 seconds: adjust intensity until a strong vibrating sensation is felt. Than adjust intensity level (during contraction phase) until a smooth tetanic contraction occurs. Discontinue treatment as soon as contraction strength drops.

#### 7.4. SEQUENTIAL PROGRAMME (TONIC FIBRES) (MF,TENS)

1. General: Complete programme for specific muscle strengthening, composed of a 'warming-up', 2 training phases, multiple intervals and 'cooling-down'.
2. Indication: In conditions, where an increment of the strength of non-injured and/or fully recovered tonic (slow-twitch) muscles is required.
3. Treatment time: 22 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 2p, 2,5 kHz, AMF 10 Hz, 1,5 min, followed by an interval of 1 min.  
BF.SYM., 400  $\mu$ s, 20 Hz, contraction rhythm of 1/6/1/8, 4 min., followed by an interval of 2 min.  
4 sequences of each 2 min. BF.SYM., 400  $\mu$ s, 10 Hz, frequency sweep 25 Hz, sweep cycle 12/12, followed by interval of 1 min.  
IF, 2,5 kHz, AMF10 Hz, 1,5 min.
6. Sensation: Sequences of relatively short, not gradually increasing/decreasing stimulation, followed by sequences with alternating rising and falling stimulation. Interrupted by intervals of different duration. Strong vibrations during first and final part.
7. Intensity level: First 90 seconds: adjust intensity until a strong vibrating sensation is felt. Than adjust intensity level (during contraction phase) until a smooth tetanic contraction occurs. Discontinue treatment as soon as contraction strength drops.

#### 7.5. SINGULAR PROGRAMME (MF)

1. General: Simple programme for non-specific muscle strengthening.
2. Indication: In conditions where an increment of muscle force is required for non-injured and/or fully recovered muscles.
3. Treatment time: 15 minutes
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: 10 sequences of each 1 min. IF, 2,5 kHz, AMF 20 Hz, frequency sweep 65 Hz, sweep cycle 1/1, contraction rhythm 1/4/1/14, followed by an interval of 30 seconds.
6. Sensation: Sequence of gradually increasing and decreasing stimulation (fast ramp-up and fast ramp-down) followed by a long interval time.
7. Intensity level: Adjust intensity level (during contraction phase) until a smooth tetanic contraction occurs. Do not increase intensity any further during operation.

#### 7.6. PROGRAMME BASED ON 'RUSSIAN STIMULATION'

1. General: Intensive programme for muscle strengthening by means of interrupted medium frequency current. The training itself is started with and terminated by a muscle relaxation programme ('warming-up' and 'cooling-down').
2. Indication: In conditions where an increment of muscle force is required for healthy or fully recovered muscles.
3. Treatment time: 23 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 2,5 kHz, AMF 10 Hz, 1,5 min.  
8 sequences of each 1½ min. RS, 2,5 kHz, 50 Hz, ratio 1:1, contraction rhythm 1/7/1/9, followed by an interval of 1 minute.  
IF, 2,5 kHz, AMF 10 Hz, 1,5 min.
6. Sensation: Sequences of fast, powerful contractions, interrupted by intervals. Strong vibrations at the start and at the end of the programme.
7. Intensity level: First 90 seconds: adjust intensity until a strong vibrating sensation is felt. Than adjust intensity level (during contraction phase) until patient's tolerance level. Discontinue treatment if contraction strength weakens.

**8.1. BURSITIS (MF)**

1. General: Programme for the improvement of the blood flow by means of medium frequency current (bipolar interferential), as part of a treatment protocol to reduce the inflammatory state.
2. Indication: Subacute and chronic bursitis, e.g. bursitis prae-patellaris.
3. Treatment time: 14 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 2p, 2,5 kHz, AMF 80 Hz, frequency sweep 50 Hz, sweep cycle 6/6, 14 min.
6. Sensation: Strong, continuous stimulation.
7. Intensity level: Clearly perceptible up to strongly perceptible.

**8.2. EPICONDYLITIS (TENS)**

1. General: Programme for the improvement of the blood flow by means of TENS current, as part of a treatment protocol to reduce the symptoms of a 'tennis elbow'.
2. Indication: Subacute and chronic tendinitis and/or periostitis of the epicondylus, e.g. epicondylitis lateralis humeri.
3. Treatment time: 15 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 75  $\mu$ s, 80 Hz, frequency sweep 20 Hz, sweep cycle 6/6, 15 min.
6. Sensation: Comfortable, fast rhythmic stimulation.
7. Intensity level: Strongly perceptible up to patient's tolerance level.

**8.3. TENDINITIS (MF)**

1. General: Programme for the improvement of the blood flow by means of a medium frequency current (bipolar interferential), as part of a treatment protocol to reduce the inflammatory state of a tendon.
2. Indication: Subacute and chronic tendinitis.
3. Treatment time: 14 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 2p, 2,5 Hz, AMF 80 Hz, frequency sweep 80 Hz, sweep cycle 1/30, 14 min.
6. Sensation: Intruding, strong vibrations.
7. Intensity level: Clearly perceptible up to strongly perceptible.

**8.4. CAPSULITIS (MF)**

1. General: Programme for the improvement of the blood flow by means of medium frequency currents (four polar interferential), as part of a treatment protocol to reduce the inflammatory state of the articular capsule.
2. Indication: Subacute and chronic capsulitis, e.g. capsulitis of the shoulder joint.
3. Treatment time: 15 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 4p, 4 kHz, AMF 100 Hz, frequency sweep 80 Hz, sweep cycle 6/6, 15 min.
6. Sensation: Tolerable, fast rhythmic stimulation.
7. Intensity level: Clearly perceptible up to strongly perceptible.

**9.1. STRESS INCONTINENCE (TENS)**

1. General: Programme for improving the strength of the striated pelvic floor muscles by means of a comfortable TENS current.
2. Indication: Stress incontinence.
3. Treatment time: 15 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 150  $\mu$ s, 35 Hz, contraction rhythm 0/7/0/25, 15 min.
6. Sensation: Comfortable, but clearly perceptible contraction of the pelvic floor, followed by relaxation.
7. Intensity level: Adjust intensity level (during contraction phase) until a smooth, clearly perceptible contraction of the pelvic floor occurs. The intensity is not further increased during operation.

**9.2. URGE INCONTINENCE (TENS)**

1. General: Programme for the inhibition of the bladder contraction and/or normalization of the tone by means of a comfortable TENS current.
2. Indication: Urge incontinence.
3. Treatment time: 15 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 100  $\mu$ s, 5 Hz, 15 min.
6. Sensation: Comfortable tingling in a low beat frequency.
7. Intensity level: Adjust the intensity level until clearly perceptible (sensoric stimulation).

**9.3. MIXED INCONTINENCE (TENS)**

1. General: Programme, in which successively pelvic floor muscles are activated and the bladder contraction is inhibited.
2. Indication: Mixed incontinence.
3. Treatment time: 20 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 150  $\mu$ s, 35 Hz, contraction rhythm 0/7/0/25, 10 min.  
BF.SYM., 100  $\mu$ s, 5 Hz, contraction rhythm 0/7/0/25, 10 min.
6. Sensation: Clearly perceptible contraction followed by a comfortable tingling sensation.
7. Intensity level: First adjust the intensity level a clearly perceptible contraction of the pelvic floor occurs (motoric level). The intensity during the second phase should be set at clearly perceptible (sensoric stimulation).

### 10.1. TROPHIC IMPROVEMENT (TENS)

1. General: Programme to improve the trophic of the tissues in the shoulder.
2. Indication: Subluxation (and pain) of the shoulder due to a C.V.A.
3. Treatment time: 90 minutes (the training must be built up from 1,5 to 6 hours a day).
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 250  $\mu$ s, 35 Hz, contraction rhythm 1/10/0/12, 90 min.
6. Sensation: Increasing stimulation followed by a long interval.
7. Intensity level: Adjust intensity level until an elevation in combination with an abduction and extension of the Humerus occurs.

### 10.2. CONTRACTURE RESTRICTION (TENS)

1. General: Programme for the restriction of contractures by means of stimulation of the antagonist.
2. Indication: Flexion contracture of the forearm due to a C.V.A.
3. Treatment time: 30 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 250  $\mu$ s, 35 Hz, contraction rhythm 1/4/1/0, 30 min.
6. Sensation: Short increasing and decreasing stimulation.
7. Intensity level: Adjust intensity level until a clearly perceptible contraction occurs which must be pain free.

### 10.3. TONE NORMALIZATION HYPERTONIA (TENS)

1. General: Programme for tone normalization by means of stimulation of the extensors of the forearm and hand for function improvement.
2. Indication: Hypertonia in the flexors of the forearm due to a C.V.A.
3. Treatment time: 10 minutes (the training must be built up from 10 to 20 min.).
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 250  $\mu$ s, 15 Hz, contraction rhythm 1/8/1/0, 10 min.
6. Sensation: Increasing and decreasing stimulation.
7. Intensity level: Adjust intensity level until a clearly perceptible contraction occurs which must be pain free.

**11.1. UPPER EXTREMITY AND TRUNK (TENS)**

1. General: Trainingsprogramme after a long period of immobilisation and in-activity, composed of a 'warming-up', training phase and 'cooling down'.
2. Indication: Muscle rehabilitation with the purpose to recover the volume of the muscle(group) by means of tetanic contractions.
3. Treatment time: 45 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 200  $\mu$ s, 8 Hz, 10 min.  
BF.SYM., 200  $\mu$ s, 45 Hz, contraction rhythm 1/4/0/15, 30 min.  
BF.SYM., 200  $\mu$ s, 4 Hz, 5 min.
6. Sensation: Sequence of gradually increasing and decreasing stimulation.
7. Intensity level: Adjust the intensity level until the desired contraction occurs. Discontinue treatment if contraction strength weakens.

**11.2. LOWER EXTREMITY (TENS)**

1. General: Trainingsprogramme after a long period of immobilisation and in-activity, composed of a 'warming-up', training phase and 'cooling down'.
2. Indication: Muscle rehabilitation with the purpose to recover the volume of the muscle(group) by means of tetanic contractions.
3. Treatment time: 45 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 350  $\mu$ s, 8 Hz, 10 min.  
BF.SYM., 350  $\mu$ s, 45 Hz, contraction rhythm 1/4/0/15, 30 min.  
BF.SYM., 350  $\mu$ s, 4 Hz, 5 min.
6. Sensation: Sequence of gradually increasing and decreasing stimulation.
7. Intensity level: Adjust the intensity level until the desired contraction occurs. Discontinue treatment if contraction strength weakens.

**13.1. ARTHRITIS**

1. Indication: Arthritis
2. Duty Cycle: 100 %
3. Treatment time: 10 minutes.
4. Intensity: 1.0 W/cm<sup>2</sup>

**13.2. BURSITIS/CAPSULITIS**

1. Indication: Bursitis/Capsulitis
2. Duty Cycle: 10 %
3. Treatment time: 8 minutes.
4. Intensity: 1.0 W/cm<sup>2</sup>

**13.3. EPICONDYLITIS**

1. Indication: Epicondylitis
2. Duty Cycle: 20 %
3. Treatment time: 4 minutes.
4. Intensity: 0.5 W/cm<sup>2</sup>

**13.4. CONTUSIONS/DISTORSIONS**

1. Indication: Contusions/Distorsions
2. Duty Cycle: 20 %
3. Treatment time: 6 minutes.
4. Intensity: 1.0 W/cm<sup>2</sup>

**13.5. LIGAMENT SPRAIN**

1. Indication: Ligament sprain
2. Duty Cycle: 50 %
3. Treatment time: 6 minutes.
4. Intensity: 1.5 W/cm<sup>2</sup>

**13.6. SCAR TISSUE**

1. Indication: Scar tissue
2. Duty Cycle: 20 %
3. Treatment time: 8 minutes.
4. Intensity: 0.8 W/cm<sup>2</sup>

**13.7. MUSCLE STRAIN**

1. Indication: Muscle strain
2. Duty Cycle: 80 %
3. Treatment time: 5 minutes.
4. Intensity: 0.4 W/cm<sup>2</sup>

**13.8. SPONDYLOSIS DEFORMANS**

1. Indication: Spondylosis deformans
2. Duty Cycle: 20 %
3. Treatment time: 10 minutes.
4. Intensity: 1.0 W/cm<sup>2</sup>

### 13.9. TENDINITIS

1. Indication: Tendinitis
2. Duty Cycle: 50 %
3. Treatment time: 5 minutes.
4. Intensity: 0.4 W/cm<sup>2</sup>

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